

## Extraction of Lithium by-products from electrolyte salt

#### using acetic acid without emitting toxic gases

Ejae Ahn, Ahyeon Kim, Leesue Ham, Cheolsoo Jung,†
Department of Chemical Engineering, University of Seoul



#### Introduction

- Recycling process of the used-battery is mainly focused on obtaining valuable metals in cathode by using inorganic acid which is harmful to the environment.
- Research on used-electrolyte recycling is insufficient, also few research has been conducted on handling the harmful gas generated through the recycling.
- In our study, we retrieved Li, F from LiPF<sub>6</sub>, HF to CH<sub>3</sub>COOLi, NaF by using acetic acid which is environmental-friendly organic acid.
- Through this process, no harmful gas is emitted and most of Li can be collected.

# Reflux condenser (2) (2) (2) (3) (4) (4) (5) (6) (7) (1) (1) (1) (2) (2) (3) (4) (4) (4) (CH<sub>3</sub>COOH(I) (1) (1) (1) (2) (2) (3) (4) (4) (CH<sub>3</sub>COOH(I) (B) (CH<sub>3</sub>COOH(I) (CH

#### **▶** Reaction Mechanism

- 1)  $LiPF_6(aq.) + CH_3COOH(l) \rightarrow HF(g) + CH_3COOLi(aq.) + PF_5(g)$
- 2)  $HF(g) + NaOH(aq.) \rightarrow NaF(s) + H_2O(l)$

 $CH_3COOH(l) + NaOH(aq.) \rightarrow CH_3COONa(aq.) + H_2O(l)$ 

 $PF_5(g) + 4H_2O(l) \rightarrow H_3PO_4(aq.) + 5HF(g)$ 

 $H_3PO_4(aq.) + 2NaOH(aq.) \rightarrow Na_2HPO_4(aq.) + 2H_2O(l)$ 


#### **▶** Reaction Condition


Reaction time: 3 hour

Reaction temperature: 80 °C, 100 °C

Reaction molar ratio:  $LiPF_6: CH_3COOH = 1:2, 1:4$ 

# Results & Discussion > Flask (1)






Wavenumber (cm<sup>-1</sup>)


FT-IR spectra of flask 1

# Table.1 FT-IR data of flask 1

Table.1 FI-IN data of Hask 1				
Acetic acid $(cm^{-1})$	Experimental (cm <sup>-1</sup> )			
1759	1564			
1709	1460			

- Upon addition of Li salt, a new peak is expected to appear at lower wavenumbers because of the  $C = 0 \cdots Li^+$  interaction.<sup>[1]</sup>
- :  $LiPF_6(aq.) + CH_3COOH(l)$  $\rightarrow HF(g) + CH_3COOLi(aq.) + PF_5(g)$
- Even if the reaction conditions change, the same substances can be seen from the experiments.





#### ► TGA Analysis

Reference CH3COOLi

#### [Theory]

- 200~300 °C : 2CH<sub>3</sub>COOLi · nH<sub>2</sub>O (100%)  $\rightarrow$  2CH<sub>3</sub>COOLi (100 ×  $\frac{66}{66+18n}$  %) + 2nH<sub>2</sub>O (100 ×  $\frac{18n}{66+18n}$  %)
- 300~500 °C: 2CH<sub>3</sub>COOLi  $\left(100 \times \frac{66}{66+18n}\%\right) \rightarrow \text{Li}_2\text{CO}_3 \left(100 \times \frac{66}{66+18n} \times \frac{74}{132}\%\right) + \text{C}_3\text{H}_6\text{O} \left(100 \times \frac{66}{66+18n} \times \frac{58}{132}\%\right)$

[80°C, LiPF<sub>6</sub>: acetic acid = 1:4]

- $CH_3COOLi \cdot 1.14H_2O (90\%) \rightarrow CH_3COOLi (73\%) + 1.14H_2O (17\%)$
- 2CH<sub>3</sub>COOLi (73%)  $\rightarrow$  Li<sub>2</sub>CO<sub>3</sub> (49%) + C<sub>3</sub>H<sub>6</sub>O (24%)


[80°C, LiPF<sub>6</sub>: acetic acid = 1:2]


- $CH_3COOLi \cdot 1.22H_2O (94\%) \rightarrow CH_3COOLi (75\%) + 0.92H_2O (19\%)$
- $2CH_3COOLi (75\%) \rightarrow Li_2CO_3 (40\%) + C_3H_6O (33\%)$
- => The number of  $H_2O$  attached to  $CH_3COOLi$  changed with the reaction ratio of  $LiPF_6$  and  $CH_3COOH$ .
- Lithium acetate is very useful in DNA analysis and is used as a buffer for gel electrophoresis of DNA and RNA.
- Lithium carbonate is an important industrial material. Its main use is a precursor for compounds
  used in lithium-ion batteries. Cathode and electrolyte are made with lithium carbonate

# Conclusions

- This study examined an environmental-friendly treatment method of lithium and harmful gas in electrolytes recycling process.
- Through this process, Li of  $LiPF_6$  was obtained in the form of  $CH_3COOLi$  and F of HF was obtained in the form of NaF.
- Li<sub>2</sub>CO<sub>3</sub> and Na<sub>2</sub>O can be obtained by thermal processing of CH<sub>3</sub>COOLi and CH<sub>3</sub>COONa.

## **Flask (2)**






#### ► FT-IR Analysis

Table.2 FT-IR data of flask 2

	Prediction	Reference(cm <sup>-1</sup> )	Experimental(cm <sup>-1</sup> )	Verification	Note
	CH <sub>3</sub> COONa	1561	1560	0	$C = O \cdots Na$
		1410	1406	Ο	interaction
	NaF	No strong peak	1	_	TGA analysis
	Na <sub>2</sub> HPO <sub>4</sub>	856	924	X	
		942	1013	X	$HPO_4^{2-}$ [3]
		1057	1045	X	

- NaF existence can't be checked by FT-IR -> Needs another method to verify its structure.
- Need a detailed plan for reaction condition of  $LiPF_6$ :  $CH_3COOH = 1:1 -> Check PF_5$  behavior.



#### ► TGA Analysis

[80°C, LiPF<sub>6</sub>: acetic acid = 1:4]

- $2CH_3COONa (40\%) \rightarrow Na_2CO_3 (26\%) + (CH_3)_2CO (14\%)$
- $Na_2CO_3$  (26%)  $\rightarrow Na_2O$  (15%) +  $CO_2$  (11%)
- Approximately 20% of weight loss can't be defined.

[100°C, LiPF<sub>6</sub>: acetic acid = 1:4]

- $CH_3COONa \cdot H_2O (56\%) \rightarrow CH_3COONa (48\%) + H_2O (8\%)$
- $2CH_3COONa (48\%) \rightarrow Na_2CO_3 (33\%) + (CH_3)_2CO (15\%)$
- $Na_2CO_3$  (33%)  $\rightarrow Na_2O$  (20%) +  $CO_2$  (13%)
- ∴ Approximately 40% NaF, theory ≈ experiment
- It is inevitable that  $PF_5$  is formed through the reaction between  $LiPF_6$  and  $CH_3COOH$ , but we couldn't find out the trace of  $PF_5$  through FT-IR analysis.
- Utilization of by-products
  - Sodium fluoride for cavity prevention and toothpaste.
  - Sodium acetate for carbon source for culturing bacteria by ethanol precipitation.

#### References

[1] Kento KimuraJoh, MotomatsuYoichi, TominagaYoichi Tominaga. "Correlation Between Solvation Structure and Ion-Conductive Behavior of Concentrated Poly(ethylene carbonate)-Based Electrolytes", The Journal of Physical Chemistry C

- 120(23) [2] 호영소 비수 전해액 및 이를 포한하는 리튬 이차 전지 10-2014-0098925 2013/01/31
- [2] 홍연숙. 비수 전해액 및 이를 포함하는 리튬 이차 전지. 10-2014-0098925,2013/01/31 [3] Saba Daneshgar, Arianna Callegari(2018), "Impact of pH and Ionic Molar Ratios on Phosphorous Forms Precipitation and Recovery from Different Wastewater Sludges", The journal of Resources.